



# Tray- and drilling protocol guide – Overlay 1






## ASTRA TECH Implant System™ EV drilling protocol overview

There is one drilling protocol recommended for preparation of the spongy bone with two options for the cortical bone.





### Step Drill

Used for site preparation – allows for a stepped osteotomy that provides guidance during the drilling process.  
Length: available in 6–13 mm and 6–17 mm.




### Cortical Drill – straight implants

Used for preparation of the cortical layer.  
**Note:** There are special cortical drills for implants with 6 mm length.



### Optional X/V osteotomy preparation



**Note:** For conical implants, color refers to the diameter of the implant body.