

ATLANTIS™ ISUS implant compatibility chart



Manufacturer/ Company	Implant system	Implant level	Abutment level	ATLANTIS™ ISUS prosthetic screw availability	Screw Order REF number from Manufacturer/ Company**
BioComp	INTERNAL				
	D4 Implant Level	x			BC-1819-225220
	D5.2 Implant Level	x			BC-1819-225220
	EXTERNAL				
	D4/D5.2 Abutment level		x		BC-1824-282015
BioHorizons	INTERNAL				
	Internal/Tapered/Single 3.5, 4.5, 5.7	x		x	PXAS
	Internal D3.5/D4.5/D5.7 Abutment for screw		x	x	PXAS
	D3.5/D4.5/D5.7 Angled Abutment for screw		x	x	PXAS
	EXTERNAL				
	External Implant 3.5	x			130-300
	External Implant 4.0, 5.0, 6.0	x			140-300
	External D3.5/D4/D5 Abutment for screw		x		222-110
Biomet 3i	INTERNAL				
	Certain™ 3.25, 4.0, 4/3, 5.0, 5/4, 6.0	x		x	IUNIHT
	Certain™ Prevail™ 3/4/3 (3.4), 4/5/4, 5/6/5, 4/3, 5/4, 6/5	x		x	IUNIHT
	Certain™ XP 4/5, 5/6	x		x	IUNIHT
	EXTERNAL				
	MicroMiniplant™ 3.25	x		x	UNIHT
	Miniplant™ 3.25	x		x	UNIHT
	Osseotite™ 3.75, 4.0, 5.0, 6.0	x		x	UNIHT
	Osseotite™ XP 4/5, 5/6	x		x	UNIHT
	Standard Abutment		x		GSH30
	Conical Abutment D3.4		x		GSH30
	Conical Abutment D4.1/D5		x		GSH30
	TG Hex Abutment		x		GSH30
	IOL Abutment		x		GSH30
Low Profile Abutment		x	x	LPCTSH	
Camlog	INTERNAL				
	Screw-line Implant 3.3	x		x	J4005.1601
	Screw-line/Root-line Implant 3.8, 4.3	x		x	J4005.1601
	Screw-line/Root-line Implant 5.0, 6.0	x			J4005.2001
	D3.3/D3.8/D4.3 Abutment level		x		J4005.1602
	D5/D6 Abutment level		x		J4005.2002
	Vario SR Abutment D3.8/D4.3, D5/D6		x		J4005.2004
	CONOLOG D3.3, D3.8, D4.3 Implant Level	x			C4005.1601
	CONOLOG D5 Implant Level	x			C4005.2001
DENTSPLY Implants	INTERNAL				
	ANKYLOS® C/X 3.5, 4.5, 5.5, 7.0	N/A			N/A
	ANKYLOS® Balance Base Narrow D4.2		x	x	3105 6022
	ANKYLOS® Balance Base D5.5		x	x	3105 6022
	OsseoSpeed™/OsseoSpeed™ TX 3.0 – Yellow*	x*			24189
	OsseoSpeed™/OsseoSpeed™ TX 3.5, 4.0 – Aqua*	x*		x	24449
	OsseoSpeed™/OsseoSpeed™ TX 4.5, 5.0 – Lilac*	x*		x	24209
	OsseoSpeed™ 3.5 – (Previous)	x*			24208
	OsseoSpeed™ 4.0 – Blue (Previous)	x*			22486
	ATIS Uni Abutment EV		x	x	25481
	ATIS 20° UniAbutment 3.5/4.0, 4.5/5.0		x	x	22435
	ATIS 45° UniAbutment 3.5/4.0, 4.5/5.0		x	x	22435
	ATIS Angled Abutment 20°		x	x	22435
	XiVE® S 3.0	x		x	46-4304
	XiVE® S 3.4, 3.8, 4.5, 5.5	x		x	46-4305
	XiVE® MP/TG 3.8, 4.5, 5.5		x	x	46-4300
	Frialit™ 3.4, 3.8, 4.5, 5.5	x		x	46-4305
	IMZ TwinPlus D3.3 Implant Level (After 1995)	x			51-1209
	IMZ TwinPlus D4 Implant Level (After 1995)	x			51-1208
	EXTERNAL				
	IMZ TwinPlus D3.3 Abutment Level		x		51-1209
	IMZ TwinPlus D4 Abutment Level		x		51-1208

* Implant-level placement is not recommended on ASTRA TECH Implant System.

** Note: We recommend to check screw REF number with the manufacturers, excluded DENTSPLY Implants.

Note: ATLANTIS ISUS assortment can vary from country to country. Please check ATLANTIS WebOrder for availability. All trademarks, company names and implant designs are the property of their respective owners.

Manufacturer/ Company	Implant system	Implant level	Abutment level	ATLANTIS™ ISUS prosthetic screw availability	Screw Order REF number from Manufacturer/ Company**
Dyna	INTERNAL				
	D3	x			82ESO
	D3.6	x			31AS1
	Octa Lock D4	x		x	81AS0
	Octa Extension Level		x		82ESO
Keystone Dental	INTERNAL				
	PrimaConnex™ SD 3.3, 3.5	x			45060K
	PrimaConnex™ RD 4.0, 4.1	x			45060K
	PrimaConnex™ WD 5.0	x			45060K
	Genesis 3.8, 4.5, 5.5, 6.5	x			45060K
	Renova SDI	x			43116
	Renova RDI	x			43116
	STAGE 1RDS Abutment level		x		S2408
	EXTERNAL				
	Restore™ SD 3.3, 3.4	x			R9203-34-48K
	Restore™ RD 3.75, 4.0	x			R9202-40-48K
	Restore™ WD 5.0, 5.5, 6.0	x			R9203-50-48K
	MultiUnit RD		x		45164K
	Detachable Abutment		x		R9460-48K
Nobel Biocare	INTERNAL				
	NobelActive™ NP 3.5	x		x	36917
	NobelActive™ RP 4.3, 5.0	x		x	36918
	NobelReplace™ NP 3.5	x		x	36818
	NobelReplace™ RP 4.3	x		x	29475
	NobelReplace™ WP 5.0	x		x	29475
	NobelReplace™ 6.0	x		x	29475
	EXTERNAL				
	Brånemark System™ NP 3.3	x			29282
	Brånemark System™ RP 3.75, 4.0	x		x	29283
	Brånemark System™ WP 5.0, 5.5	x			29284
	Multi-Unit Abutment RP	x		x	29285
	Multi-Unit Abutment WP		x		29286
	Standard Abutment		x		32983
EsthetiCone Abutment		x		32981	
Steri-Oss PME Abutment		x		2346	
Straumann	INTERNAL				
	Bone Level 3.3 NC (Narrow CrossFit™)	x		x	025.2900
	Bone Level 4.1, 4.8 RC (Regular CrossFit™)	x		x	025.4902
	NCC	x			048.313
	Standard/Standard Plus 4.8 RN (Regular Neck)	x		x	048.356
	Standard/Standard Plus 4.8 WN (Wide Neck)	x		x	048.356
	Standard Plus 3.3 NN (Narrow Neck)	x			049.177
	EXTERNAL				
	Monotype 8° Conus		x	x	048.350
	RN Abutment Level		x	x	048.350
	RN Abutment Level (Old)		x		048.350
	WN Abutment Level		x	x	048.350
	Bone Level Multi-Base Abutment D3.5, D4.5		x	x	025.0900-04
	Bone Level Multi-Base Abutment D6.5		x		025.0900-04
Bone Level Multi-Base Angled Abutment		x	x	025.0900-04	
Zimmer Dental	INTERNAL				
	Tapered Screw-Vent™ 3.5, 4.5, 5.7	x		x	MHLAS
	Screw Vent™ 3.3, 3.7, 4.5	x		x	MHLAS
	SwissPlus™ 3.8 Narrow	x			MHLAS
	SwissPlus™ 4.8	x			GPCAS
	Advent D4.5	x			AVHLS
	EXTERNAL				
	Spline 3.25, 4, 5mm	x			1537
	Screw Vent™ TSI Abutment		x		TSF
	Tapered Abutment		x	x	SCTS
	Spectra Cone Abutment		x	x	SCTS
Spline Shouldered Abutment D4.5, 6.5		x		1621	

** Note: We recommend to check screw REF number with the manufacturers, excluded DENTSPLY Implants.
 Note: ATLANTIS ISUS assortment can vary from country to country. Please check ATLANTIS WebOrder for availability.
 All trademarks, company names and implant designs are the property of their respective owners.



ATLANTIS™ ISUS turnaround times

ATLANTIS ISUS implant suprastructures are available for all major implant systems, in titanium and cobalt-chrome. Turnaround times vary depending on material.

When the case is received at the DENTSPLY Implants production facility and the manufacturing process has started, your ATLANTIS WebOrder account will be updated with the case estimated shipment date. Standard turnaround times are as follows:

ATLANTIS™ ISUS in-house production time (days)		
	Titanium	Cobalt-chrome
ATLANTIS™ ISUS Bridge	5	5
ATLANTIS™ ISUS Hybrid	5	5
ATLANTIS™ ISUS Bar	5	5
ATLANTIS™ ISUS 2 in1	10	N/A

Determine return date

1. Select your choice of delivery method and calculate the number of days it will take to ship the case from the dental laboratory to the DENTSPLY Implants production facility.
2. Add the days the case will be in-house at the production facility. Calculation of in-house production time starts when the case is received at the production facility before noon.
3. DENTSPLY Implants will ship the case free of charge via next day delivery.

Note:

- Design approvals to process case must be granted by 6 p.m. CET.
- Material, on-hold and case information issues must be resolved by midnight CET.
- Weekends and holidays are excluded.
- Time needed for **Design approval** or resolving **On hold** information may delay the case.

Example

Order and all case materials for an ATLANTIS ISUS Bridge in titanium is received on Monday before noon. Design for approval will be sent out.

If approval takes place the same day by 6 p.m. CET, the case will be shipped to the dental laboratory on Friday.

Additional information

The **ATLANTIS ISUS prosthetic screw** is offered for selected major implant systems. The screws can be ordered together with the ATLANTIS ISUS implant suprastructure.

The ATLANTIS ISUS implant suprastructures are always compatible with the implant/abutment screw from respective manufacturers system. If the ATLANTIS ISUS prosthetic screw is not available, the screws can be ordered from the implant system manufacturer. As a service all screws are listed with the manufacturers REF number in the ATLANTIS ISUS compatibility chart.

Note: We recommend to check the REF number with the implant system manufacturer.

For the most current version of the ATLANTIS ISUS implant compatibility chart and turnaround times, visit www.dentsplyimplants.com or your ATLANTIS WebOrder account.

